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---

**Instructions for use**  
**ACUVUE® contact lenses for daily use (frequent replacement)**

Brand	Material	Packing Solution	Usage Plan	ANVISA license
ACUVUE® 2	etafilcon A	1	14 days (for daily use) 6 nights/7 days (for extended use)	80145909100
ACUVUE® OASYS with HYDRACLEAR® PLUS 3	senofilcon A	2	14 days (for daily use) 6 nights/7 days (for extended use)	80145909101
ACUVUE® VITA 3	senofilcon C	2	1 month (for daily use)	80145909103
ACUVUE® VITA for ASTIGMATISM 3	senofilcon C	2	1 month (for daily use)	80145909104
ACUVUE® OASYS for ASTIGMATISM with HYDRACLEAR® PLUS 3	senofilcon A	2	14 days (for daily use) 6 nights/7 days (for extended use)	80145909102
ACUVUE® OASYS MULTIFOCAL 3	senofilcon A	2	14 days (for daily use) 6 nights/7 days (for extended use)	80145909105
ACUVUE® OASYS with Transitions™ 3 4	senofilcon A	2	14 days (for daily use)	80145909106

**Manufactured by:** Johnson & Johnson Vision Care, Inc. - 7500 Centurion Parkway - Jacksonville, Florida 32256 - USA

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**Resp. Tech.:** Carla Carvalho de Silveira CRQ No. 04343300

**Validity (year / month) and Lot:** see package side

**Sterile:** steam sterilized product.

**Content:** contact lenses

SALE UNDER REFRACTIONAL MEDICAL PRESCRIPTION.

REPROCESSING PROHIBITED.

Store the product in its original packaging at room temperature and away from direct sunlight.

Do not use if packaging is violated.

Discard after recommended use.

See Instructions for Use for Warnings, Precautions and other information.

Caption: **Packaging** : 3 Solution : 1 Borate Saline Buffered Saline 2 Solution Borate Buffered with Cellulose Methyl Ether **Material** : 3 The lens material contains silicone and meets Class 1 UV absorption standards with less than 1% UVB radiation transmissibility (280 -315 nm) and 10% UVA radiation (315-380 nm). All other ACUVUE® products meet Class 2 UV absorption standards with transmissibility of less than 5% UVB radiation and 50% UVA radiation 4. The lens material contains a photochromic additive that dynamically absorbs visible light over a range of 380 nm to 780 nm upto a minimum transmissibility of 84% in the inactivated (closed) state to a minimum of 23% in the activated (open) state in a function of lens thickness and level of absorbed UV and high energy visible light (HEV) radiation.

<sup>1</sup> Photochromic is defined as exhibiting a reversible change in color or tone when exposed to specific types of light of sufficient intensity. In the absence of activating light, the color returns to normal color.

**Consumer Services**

0800 762 5424

[www.acuvue.com.br](http://www.acuvue.com.br)

## 1. PRODUCT DESCRIPTION and INTENDED USE

This leaflet concerns the following ACUVUE® contact lenses that are intended for daily wear (Frequent replacement) (less than 24 hours while you are awake) or extended wear (more than 24 hours, including while you sleep) which can be worn continuously for up to 7 days/6 nights and must be discarded after removal. If used in this way, no cleaning or disinfection is required. Your ophthalmologist will determine your usage plan and replacement plan.

- ACUVUE® Spherical Contact Lenses are intended for daily use and for the optical correction of nearsightedness (difficulty seeing at a distance) and hyperopia (difficulty seeing up close) in people with healthy eyes that may be up to 1.00D (diopter). ) or less of astigmatism.

- **ACUVUE® OASYS with Transitions™** contact lenses are indicated for daily optical refractive correction wear for patients who have difficulty seeing at a distance (myopia) and up close (hyperopia) who may have 1.00D or less of astigmatism and who do not have any eye disease. They are also suitable for attenuating light, as it contains a photochromic additive that dynamically absorbs visible light.

- **ACUVUE® Contact Lenses for ASTIGMATISM** are indicated for daily use for vision correction in people with healthy eyes who are nearsighted or farsighted and may have 0.50 D to 3.00 D of astigmatism.

- **ACUVUE® OASYS Contact Lenses with HYDRACLEAR® PLUS** are indicated for the optical correction of refractive ametropia (myopia and hyperopia) in aphakic and non-aphakic individuals with healthy eyes who have 1.00D or less of astigmatism. They are also indicated for therapeutic use as a dressing lens (bandage).

- **ACUVUE® OASYS Contact Lenses with HYDRACLEAR® PLUS** are also indicated for therapeutic use as a treatment lens for certain eye problems.

- **ACUVUE® OASYS MULTIFOCAL** contact lenses are indicated for daily use for vision correction, specifically, optical correction of ametropia and presbyopia, and also for therapeutic use as a dressing lens (bandage).

Your eye care professional will determine your usage schedule (how long you should wear your lenses per day) and your replacement schedule (when you should discard your lenses and wear new ones). When prescribed for frequent/planned replacement use, you can clean and disinfect lenses using only a chemical disinfection system.

Your lenses are intended for daily wear for up to 14 days (2 weeks) and extended wear for up to 6 nights / 7 days of continuous wear. It is recommended that you first be evaluated on a daily usage schedule. If successful, a gradual introduction of long-term use can be followed, as determined by your ophthalmologist.

When indicated for therapeutic use may aid the treatment of certain ocular conditions and acute conditions:

- For corneal protection from eyelid and corneal anomalies such as entropion, trichiasis, tarsal scarring ( lower eyelids ) and recurrent corneal erosion. Additionally, they are also indicated for protection whenever sutures or malformation of the ocular structure, degeneration or paralysis may result in the need to protect the cornea from repeated exposure or irritation.

- For corneal pain relief in conditions such as bullous keratopathy, epithelial erosion and abrasion, filamentary keratitis, and post- keratoplasty .

- For use as a barrier during the healing process of epithelial defects, such as chronic epithelial defects, corneal ulcer, contact dermatitis neurotrophic and neuroparalytic and chemical burns.

- For post-surgical conditions for which the use of treatment lenses is indicated, such as post-refractive surgery, lamellar grafts, corneal folds, and additional ocular surgical conditions.

- For stability and structural protection in piggyback lens implants where the cornea and associated surfaces are too irregular to allow placement of rigid gas permeable contact lenses (RGP) on the cornea. Additionally, lens wear can prevent irritation and abrasions in conditions where there are elevation differences at the host/graft junction or scar tissue.

- Lenses prescribed for therapeutic use can be worn daily or for extended periods .

All ACUVUE® contact lenses have a UV radiation filter that helps protect against harmful transmission of UV radiation to the cornea and eye.

**WARNING: UV RADIATION ABSORBING CONTACT LENSES do not replace the use of other eye protection methods that absorb UV radiation, such as protective eyewear or sunglasses, as they do not fully cover the eyes and the areas around them . You should continue to wear UV radiation absorbing eye protection.**

**NOTE: Prolonged exposure to UV radiation is one of the risk factors associated with the development of cataracts. Exposure is based on various factors such as environmental conditions (altitude, geography, cloudy skies) and personal factors (duration and nature of outdoor activities). Radiation (UV) blocking contact lenses help protect against damage caused by UV radiation. However, clinical studies have not yet been performed to demonstrate that wearing contact lenses with a radiation (UV) blocker reduces the risk of developing cataracts or other eye diseases. Consult your ophthalmologist for more information.**

## 2. CONTRAINDICATIONS (in which situations you should not use)

When wearing contact lenses for vision correction, DO NOT wear these lenses when you have any of the following conditions:

- Inflammation or infection in or around the eyes or eyelids
  - Any eye disease, injury or abnormality that affects the cornea, conjunctiva or eyelids
  - Any previously diagnosed problem that makes wearing contact lenses uncomfortable
  - Severely dry eye
  - Reduced corneal sensitivity
  - Any systemic disease that can affect the eyes or worsen with contact lens use
  - Allergic reactions on the surface of the eyes or in the tissues around them that can be induced by the use of contact lenses or worsen due to their use.
  - Any active eye infection
  - If eyes get red or irritated
- Eye irritation caused by allergic reactions to ingredients in contact lens solutions (eg, cleaning and disinfecting solution, lubricating drops, etc.). These solutions may contain chemicals or preservatives (such as mercury, thimerosal, etc.) to which some people may develop an allergic response.

Your eye care professional may prescribe contact lenses for therapeutic use to aid the healing process of certain eye conditions, which may include those listed above.

### **3. WARNINGS - What You Should Know About Using Contact Lenses:**

SALE UNDER REFRACTIONAL MEDICAL PRESCRIPTION.

PROHIBITED REPROCESSING PRODUCT.

EYE PROBLEMS, INCLUDING CORNEAL ULCER, CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION. IF YOU FEEL:

#### **Instructions for use**

**ACUVUE® contact lenses for daily use (frequent replacement) and use prolonged the**

- Eye discomfort,
- Excessive tearing,
- Change in vision,
- Loss of vision,
- Redness in the eyes, or
- Any other eye problems,

#### **Instructions for use**

**ACUVUE® contact lenses for daily use (frequent replacement) and use prolonged the**

YOU MUST REMOVE YOUR LENSES IMMEDIATELY AND CONTACT YOUR OPHTHAMOLOGIST.

When prescribed by your ophthalmologist for daily wear (ie, instruct you to remove your lenses at the end of each day), you should not wear your contact lenses while you sleep (outside of the prescribed wearing schedule). Clinical studies have shown that there is a serious risk of eye problems when the lenses are worn at night.

Studies have shown that contact lens wearers and smokers are at greater risk of eye problems than non-smokers.

Problems with contact lenses or lens care products can result in serious eye damage.

The correct use of your contact lenses and lens care products, including lens cases, are essential to the safe use of this product.

The overall risk of serious eye problems can be reduced by carefully following lens care guidelines, including lens case cleaning.

#### **Specific warnings:**

##### ***Activity in water:***

**Instructions for use:** Do not expose your contact lenses to water while you are wearing them.

**Warning:** Water may contain microorganisms that can lead to severe infection, loss of vision or blindness. If your lenses are submerged in water while participating in water sports or swimming in pools, hot tubs, lakes and oceans, you must discard them and replace them with new lenses. Ask your eye care professional for recommendations on wearing lenses during any activity involving water.

##### ***Immersion and storage of your lenses:***

**Instructions for use:** Do not reuse disinfectant solutions to store your contact lenses.

**Warning:** Do not reuse or top up the old disinfectant solution contained in your case as reused solutions reduce the effectiveness of lens disinfection and can lead to severe infection, loss of vision or blindness.

##### ***Multipurpose solution discard date:***

**Instructions for use:** Discard any remaining disinfectant solution after the indicated time. Disposal time indicates the safe period in which you can use the disinfectant solution after opening the package, which may not be the same as the product's expiration date, which brings the product's expiration date before

opening the bottle.

**Warning:** Using your solution for disinfection beyond the disposal date can result in contamination of the solution and can lead to severe infection, loss of vision or blindness. To avoid contamination, DO NOT touch the tip of the container to any surface. Replace cap after use. To avoid contaminating your solution, DO NOT transfer the solution to other containers or containers.

**Friction and rinsing of lenses:**

**Instructions for use:** To properly disinfect your contact lenses, rub and rinse your contact lenses according to the time recommended in the instructions for use of your disinfectant solution.

**Warning:** Rub and rinse contact lenses for the recommended time to help prevent serious eye infections. Never use water, saline (saline solution), or moisturizing drops to disinfect lenses. These solutions will not disinfect lenses. Failure to use the recommended disinfectant can lead to serious infection, loss of vision or blindness.

**Care of the lens storage case:**

**Instructions for use:** Empty and clean your contact lens case by rubbing it using sterile disinfectant solutions or contact lens cleaning solutions. NEVER use water. Cleaning should be followed by rinsing with a sterile disinfectant solution (never use water) and it is recommended to dry the lens case with a clean paper towel. Never air dry or cover the lens case without using any additional cleaning methods. If drying in air, ensure that no residual solution remains in the case. Replace your lens case according to the instructions given by your eye care professional or on the label that came with your lens case.

Contact lens cases can be a source of bacterial growth.

**Warning:** Do not store your lenses or wash your lens case with water or any non-sterile solution. Use only appropriate multipurpose solution so there is no contamination of your lens or lens case. The use of non-sterile solution can lead to serious infection, loss of vision or blindness.

<sup>2</sup> New England Journal of Medicine, September 21, 1989; 321 (12), pp. 773-783

#### 4. PRECAUTIONS

**For the health of your eyes,** it is important to carefully follow the handling, placement, removal and use instructions contained in this guide, as well as any extra care prescribed by your ophthalmologist ( **see the sections “Lens handling and placement” and “Lens use ”** ).

**General Precautions:**

If you wear contact lenses to correct presbyopia using monovision, you may not be able to get the best corrected eyesight for near or far vision. Visual needs are different for different people, so your ophthalmologist should work with you to select the most appropriate type of lens for your needs.

- **Always** contact your eye care professional before using any medication in your eyes.
- **Be alert** as certain medications such as antihistamines, decongestants, diuretics, muscle relaxants, tranquilizers and those for motion sickness can cause dry eyes, increased lens perception (sensation of the presence of lenses in the eyes) or blurred vision. Always tell your eye care professional if you have problems with your lenses while using these medications. Depending on your symptoms, your ophthalmologist may recommend applying moisturizing drops suitable for use with soft contact lenses, or may recommend that you stop wearing contact lenses when using these medications.
- **Be alert** as oral contraceptive users may develop changes in visual quality and comfort while wearing contact lenses.

Do not change your contact lens type (eg make, model, technology) or parameters (eg diameter, base curve, diopter, etc.) without consulting your ophthalmologist.

- **Be aware** that wearing the ACUVUE® OASYS lens with Transitions™ on only one eye is not recommended as it can disrupt your ability to perceive depth and movement of objects. It can also cause an aesthetic concern.
- **Never** use ACUVUE® OASYS with Transitions™ as protection against artificial light such as lamps, lasers, etc.
- **Do not** use ACUVUE® OASYS with Transitions™ while you sleep.
- **Never** look directly into the sun or an eclipse with or without ACUVUE® OASYS contact lenses with Transitions™.
- **Be aware** that there is no data to support the safety or performance of a patient aged 50 years or older with lens opacity (eg cataract ) driving using ACUVUE® OASYS with Transitions™.
- **Always** have eyeglasses that match your current diopter available for situations where you are unable to wear your contact lenses or in situations where wearing contact lenses is not recommended.

With the use of any contact lens, regular eye care visits are necessary to ensure your eyes' safety. Consult your ophthalmologist about how often you should visit him.

**Who Should Know You Are Wearing Contact Lenses:**

- **Inform** all your doctors that you wear contact lenses.
- **Always** inform your employer that you are a contact lens wearer. Some activities may require you to wear protective eye gear or may require you not to wear contact lenses.

#### 5. ADVERSE REACTIONS (Side effects) – Possible problems and what to do

### **Possible reactions:**

Be aware that reactions can occur while you are wearing contact lenses and may or may not be associated with the following symptoms:

- Burning, pinching/stinging and/or itchy eyes
- Reduced feeling of comfort with contact lenses
- Sensation of something in the eye (foreign body, scratched area)
- Swelling or inflammation in or around the eyes
- Redness in the eyes
- Problems with the eyelids
  
- Lachrymation
- Unusual secretions in the eyes
- Reduced visual quality
- Blurred vision
- Perception of rainbows or halos around objects
- Sensitivity to light (photophobia)
- Eye dryness

When any of the above symptoms occur, a serious eye problem may be present. **You should immediately see your eye care professional** so that the problem can be identified and treated, if necessary, to avoid serious eye damage.

### **How to recognize problems and what to do**

You should perform a simple three-part self-examination at least once a day. Ask yourself:

- How do the lenses feel in my eyes?
- How do my eyes look?
- Did I notice any changes in my vision?

If you notice any problems, you should REMOVE YOUR LENS IMMEDIATELY. If the problem or discomfort stops and the lens appears to be undamaged, clean and rinse the lens with a recommended soft contact lens care solution and reinsert the lens. If, after reinserting the lens, the problem persists, discard the lens and place a new lens in the eye.

If, after inserting the new lenses, the change continues, **REMOVE THE LENSES IMMEDIATELY AND CONTACT YOUR OPHTHALMOLOGIST.**

During therapeutic use, an adverse effect may be due to the original illness or the injury may be due to the effects of contact lens wear. There is a possibility that the existing disease or condition could worsen when a soft contact lens for therapeutic use is used to treat an already sick or injured eye. To avoid serious eye damage, you should contact your Ophthalmologist IMMEDIATELY if there is an increase in symptoms during lens wear.

DO NOT use a new lens to self-treat the problem.

### **INSTRUCTIONS FOR USE**

**For the health of your eyes**, it is important to carefully follow the handling, placement, removal and use instructions contained in this guide, as well as the extra care prescribed by your ophthalmologist. If you are unable or unable to follow the recommended care procedures, you should not attempt to wear contact lenses.

**When you buy your lenses for the first time, make sure you can put them on and off (or make sure you are advised by someone else who can take the lenses out for you properly).**

#### ***Step 1: Getting Started***

It is essential that you learn and have good hygiene in the care and handling of your new lenses.

Cleanliness is the first and most important aspect of proper contact lens care. In particular, your hands should be clean, dry and free from any soap, lotion or cream residue before handling the lenses.

Before starting:

- Always wash your hands thoroughly with a mild soap, rinse them thoroughly, and dry them with a lint-free towel before touching the lenses.
- Avoid using soaps containing moisturizing cream, lotions, or oily cosmetics before handling lenses. These substances can come into contact with the lenses and interfere with their successful use. It's best to put your lenses on before you put on makeup.

#### ***Step 2: Opening the package***

Always confirm that the parameters (eg diameter (DIA), base curve ( CB ), diopter (D), etc.) printed on the cartridge and blister are in accordance with your prescription. **DO NOT USE** if there is any discrepancy.

Always store the individual blister away from direct light.

## multiple packaging

Each multipack contains individually wrapped lenses. Each lens comes in its own packaging specifically designed to keep it sterile while the packaging is sealed.

Always store unopened containers of ACUVUE® OASYS Transitions™ out of direct sunlight.

## lens packaging

**DO NOT use if sterile packaging is broken or damaged.**

To open your individual sterile package, follow the simple steps below:

1. Shake the package and verify that the lens is floating in the solution.
2. Remove the laminated seal to gain access to the lens.
3. Place a finger on the lens (do not touch it with your fingernail) and slide it out of the packaging compartment until the lens is out of it.

**NEVER use tweezers or other equipment to remove the lens from its compartment.**

Occasionally, a lens may attach to the inside surface of the laminated seal when opened, or to the plastic in the packaging itself. It will not affect the sterility of the lens. It will still be perfectly safe to use. Carefully remove the lens and inspect it following the handling instructions.

## Tips for handling lenses

- Handle your lenses with your fingers and be careful to avoid contact with your fingernails. It's good to keep your nails short and filed.

Get into the habit of always handling the same lens first to avoid mixing them up.

Once you have removed the lens from its packaging or contact lens case, check it to make sure it is moist, clean, and free of any nicks or tears. DO NOT use the lens if there is any evidence of damage.

## Step 3: Putting the lenses in the eyes

Remember to always start with the same eye. Once you have opened the lens package, removed and examined it, follow these steps to place the lens in the eye:

1. ENSURE THAT THE LENS IS NOT IN ITS POSITION INVERTED, following the procedure:

- Place the lens on the tip of your index finger and check it. The lens should have a natural, curved, round shape. If the lens edges tend to point outward, the lens is inverted.
- Another method is to lightly press your contact lens between your thumb and forefinger. Edges should turn inwards. If the lens is inverted, the edges will turn out slightly.

1. With the lens on your index finger, use your other hand to hold the upper eyelid so you don't blink.
2. Pull down your lower eyelid, using the other fingers of your hand with the lens.
3. Look up and gently place the lens on the lower white part of the eye.
4. Release your eyelid slowly and close your eye for a moment.
5. Blink several times to center the lens.
6. Use the same technique when inserting the lens into the other eye.

There are other methods for lens placement. If the above method is difficult for you, ask your ophthalmologist about an alternative method.

## Step 4: How to check if your lenses are correctly placed

Once you have successfully inserted your lenses, ask yourself:

- Do I see well?
- How do the lenses feel in my eyes?
- Do my eyes look okay?

If after putting the lenses in, your vision is blurry, check the following:

- The lenses are not centered on the eyes (see "Step 5: How to Center the Lenses" below in this instruction for use).
- If the lens is centered, remove it (see "Removing Your Lenses") and check the following points:
  - Cosmetics or oil on the lens. Clean your lenses by rubbing them with a cleaning solution, and if the cosmetic or oil is removed, put the lenses back into the eye. If your vision remains blurry, discard the lens and insert a new lens.
  - The lens is in the wrong eye.
  - The lens is inverted (also not comfortable as it should be). See "Step 3: Putting the lens in your eyes".

If you detect that your vision is still blurry after checking the above possibilities, remove both lenses and see your eye care professional.

**Note:** If a lens is noticeably uncomfortable when putting it on or becomes uncomfortable soon after insertion, remove the lens immediately and contact your eye care professional. If your eye and lens self-examination shows any other problems, IMMEDIATELY REMOVE THE LENSES AND CONTACT YOUR OPHTHALMOLOGIST.

## Step 5: How to center the lenses

A lens, positioned on the cornea (center of the eye), very rarely moves to the white part of the eye during use. However, this can occur if the insertion and removal procedures are not carried out properly. To center a lens, do one of the following:

- Close your eyelids and lightly massage the lens into place through the closed eyelids.

OR

- Gently move the lens off-center to the cornea (center of the eye) while the eye is open, using the upper or lower lid edge to move the lens into the correct position.

## USE OF LENSES

When wearing lenses, remember the following important precautions:

### **risk conditions**

- If you use aerosol (spray) products such as hairspray, when wearing lenses, keep your eyes closed until the spray dissolves in the air.
- Avoid exposing yourself to any harmful or irritating vapors, as well as smoke, when wearing lenses.
- NEVER clean your contact lenses with tap water. This water can be contaminated with various impurities that can contaminate or damage your lenses and can cause infection or injury to the eye.

### **water activity**

- Do not expose your contact lenses to water while wearing them.

### **How to lubricate/apply solutions**

- Your eye care professional can recommend an appropriate lubricating/humectant solution to moisten (lubricate) the lenses when you are wearing them.
- Do not use saliva or anything other than the recommended solution to lubricate your lenses. Do not put lenses in your mouth.

### **sharing lenses**

- Never share your lenses with another person. Sharing lenses greatly increases the chance of eye infections.

### **Adherence to prescribed replacement and usage schedule**

- Never use your lenses beyond the time indicated in the product's instructions for use or by your ophthalmologist.
- Always dispose of used lenses as directed in the instructions for use and by your ophthalmologist.

## **REMOVING YOUR LENSES**

**CAUTION:** Always make sure the lens is in the cornea (in the center of your eye) before attempting to remove it. Confirm this by covering the other eye. If vision is blurry, the lens is in the white part of the eye or not in the eye. To locate the lens, inspect the upper eye area, looking down in the mirror while lifting the upper eyelid. Then inspect the lower area by pulling the lower eyelid down.

1. Wash, rinse and dry your hands thoroughly. You should follow the method recommended by your ophthalmologist. Below is an example of a lens removal method:

### **Clamping method:**

**Step 1.** Look up, slide the lens to the bottom of the eye using your index finger.

**Step 2.** Gently pinch the lens between your thumb and forefinger.

**Step 3.** Remove the lens.

2. Follow the instructions in the next section "Caring for your contact lens".

Note: To maintain your eye health, it is important that the lens moves freely in your eye. If the lenses stick (stop moving) in the eye, apply a few drops of the moisturizing solution recommended by your ophthalmologist. Wait until the lenses are moving freely in your eyes to remove them. If they continue not to move, you should see your eye care professional immediately.

## **TAKING CARE OF YOUR LENSES**

When you first receive your contact lenses, you will receive instructions for your cleaning and disinfecting routine as recommended by your ophthalmologist. If you do not follow these instructions for use, you may develop serious eye problems, as described in the "Warnings" section.

### **Before starting:**

- Always wash your hands with a mild soap, rinse your hands thoroughly, and dry with clean paper towels before touching your contact lenses.
- Always have at hand:
  - Disinfectant solutions or multipurpose solutions recommended by your ophthalmologist. Do not change your cleaning solution without consulting your eye care professional. NEVER use cleaning solutions recommended for rigid lenses only.
  - Clean your contact lens case.

When using hydrogen peroxide systems to care for your lenses, only use the case provided for this type of cleaning solution. This kit was specially developed to neutralize the solution. Failure to use the indicated case may cause burning, burning or eye damage.

### **Step 1: Cleaning and Disinfection (Chemistry - No heating)**

1. Always clean and disinfect the same lens first to avoid mixing them up.

2. Always follow the directions for use of your cleaning solution or the instructions recommended by your ophthalmologist.

**NEVER use heat to disinfect your lenses.**

**WARNING: Lenses disinfected with chemical cleaning solutions can absorb ingredients contained in**

**the lens and may cause eye irritation.** Thorough cleaning with sterile saline solution (or other recommended solution) before putting on contact lenses can decrease the potential for irritation.

1. Rub and rinse your lenses according to the time indicated in your multi-purpose solution's instructions for use to properly disinfect your lenses.

**WARNING:**

- Rub and rinse your lenses for the recommended amount of time to avoid severe eye infections.
- Never use water, saline, or moisturizing solutions to disinfect your lenses. These solutions will not disinfect your lenses. Failure to use the recommended disinfection solution can lead to severe infection, loss of vision and blindness.

2. Put your lenses in the correct camera in its own case to store your lenses.

**Step 2: Storage**

- To store your lenses, first disinfect them and then leave them in their closed case until you are ready to use them again.

ALWAYS keep your lenses completely covered with the recommended disinfectant solution while not in use. Long periods without a disinfectant solution will make it more difficult for your lenses to return to a hydrated state. If your lenses become stiff and dehydrated, discard them and put in new lenses.

- Always store your contact lenses in lens cases and out of direct sunlight.
- If you will not be using your contact lenses right after disinfection, ask your eye care professional about storing your contact lenses.
- Always empty and clean your contact lens case by rubbing it with your finger and using cleaning solutions. After cleaning, you should rinse your own contact lens case with a disinfectant solution (never use water) and dry it with a clean paper towel. If the case is air-drying, make sure there is no disinfectant solution residue.
- Change your lenses as recommended by your ophthalmologist or in the instruction for use.
- Contact lenses can be a source of bacterial growth.


**WARNING:** Do not clean or store your contact lens case with water or any other non-sterile solution. Only use multipurpose solutions to avoid contamination of your contact lenses. Using a non-sterile solution can cause severe infections, loss of vision or blindness.

**6. NOTIFICATION OF ADVERSE REACTIONS (Side effects)**

Any incident while using ACUVUE® Contact Lenses should be reported to the manufacturer and/or its authorized representative and/or regulatory authority in your country.

**SYMBOLS**

The following symbols may appear on the product label or box:

Symbol	Definition
	Caution, see instructions for use
	Made by or in
	Manufacturing date
	Product Expiration/Valid Date
	Product manufacturing batch number
	Steam Sterilization
	Do not use if packaging is violated
	Regarding disposable material
	ATTENTION: sale under medical prescription
	Medical device in the European Community
	Indicates a unique sterile barrier system

   Filtro UV UV Blocking	UV filter
	Authorized Representative in the European Union
	CE mark and notified body identification number
<b>DIA</b>	Diameter
<b>BC</b>	base curve
<b>D</b>	Diopter ( lens degree )
<b>CYL</b>	Cylinder
<b>AXIS</b>	Axle
<b>MAX ADD</b>	Maximum addition that can be corrected
<b>LOW/ L</b>	Addition of close "low"
<b>MID/ M</b>	Addition of close to "average"
<b>HGH/ H</b>	Addition of close " high "
	Correct lens orientation
	Incorrect lens orientation (lens inside out)
	Keep away from direct sunlight