



TECNIS®

Toric II 1-Piece IOL

ABOUT TECNIS® TORIC II IOL (Continued)

The TECNIS® Toric II 1-Piece posterior chamber lens is indicated for the visual correction of aphakia and pre-existing corneal astigmatism in adult patients with or without presbyopia in whom a cataractous lens has been removed by phacoemulsification and who desire improved uncorrected distance vision, reduction in residual refractive cylinder, and increased spectacle independence for distance vision. The device is intended to be placed in the capsular bag.

INDICATIONS AND IMPORTANT SAFETY INFORMATION FOR TECNIS® Toric 1-Piece IOL and TECNIS® Toric II 1-Piece IOL

INDICATIONS: The TECNIS Toric II 1-piece posterior chamber lenses are indicated for the visual correction of aphakia and pre-existing corneal astigmatism in adult patients with or without presbyopia in whom a cataractous lens has been removed by phacoemulsification and who desire improved uncorrected distance vision, reduction in residual refractive cylinder, and increased spectacle independence for distance vision. The device is intended to be placed in the capsular bag.

IMPORTANT SAFETY INFORMATION: Routine cataract surgery risks, which are not related to the intraocular lens ("IOL"), could include minor, temporary, or permanently impact your vision. While rare, complications could include worsening of vision, bleeding, or infection. Please discuss all risks and benefits with your Eye Care Professional before lens selection and surgery. If you have a pre-existing eye condition (including glaucoma), your vision may not be good even after cataract removal, as you will not be able to benefit fully from the insertion of the IOL. Before surgery, your Eye Care Professional will check for any eye diseases. Once the IOL has been inserted, your vision may not be good enough to perform detailed 'up-close' work without glasses. While rare, in some cases certain types of retinal treatment (e.g., retinal tear repair) may become more difficult to repair. You should avoid any activity that could harm your eye(s) while you are recovering from surgery. Your Eye Care Professional will advise which activities to avoid while you are recovering from surgery. A Toric IOL corrects astigmatism when it is placed correctly in the eye. There is a chance that the Toric IOL could be placed incorrectly during the surgery or could move within the eye following the surgery. If the Toric IOL is not placed correctly, this may lead to visual distortions. A subsequent surgery may be required to adjust the position of the IOL.


SIDE EFFECTS: In 3.4% of patients, side effects with TECNIS® Toric 1-piece IOL included the need for a second surgery to reposition the lens or repair the retina. Other side effects included swelling/thickening of the retina (2.9%), and separation of the retina from the surrounding tissues (less than 1%).

ATTENTION: Reference the Directions for Use for a complete listing of Indications and Important Safety Information.

REFERENCES: 1. Freeman W. 2017 IOL report: a global market analysis for 2016 to 2022 [market report]. St. Louis, MO: Market Scope LLC; 2017.

TECNIS is a trademark of Johnson & Johnson Surgical Vision, Inc.
© Johnson & Johnson Surgical Vision, Inc. 2022 | PP2022CT5146

CLEAR OR CRISP: HOW ABOUT BOTH?



Toric IOLs for Astigmatism

LIVE THE LIFE YOU WANT TO SEE.

TECNIS®
Toric II 1-Piece IOL

Toric II

About TECNIS® Toric II IOL: Physicians considering use of the TECNIS® Toric II IOLs should refer to the Directions for Use labeling for a complete list of indications and safety information.

See reverse for continued Indications and Important Safety Information.

TRANSFORM YOUR VISION.



Astigmatism means: the shape of your eye is more oblong than circular, which prevents you from focusing properly. **Astigmatism may affect your ability to see crisp, clear vision. With the TECNIS® Toric II IOL for astigmatism, you can LIVE THE LIFE YOU WANT TO SEE.**



If you have both **astigmatism** and **cataracts**, you may have the opportunity to address them simultaneously with a toric IOL.


FIND YOUR FOCUS.

Cataract surgery is trading one thing for another. When your cataract is removed, it is replaced by a small, soft lens called an *intraocular lens* — or *IOL* for short. You won't be able to feel or see the IOL, but it will stay inside your eye to help you see clearly.


Toric IOLs can address astigmatism as part of cataract surgery.



TECNIS® Toric II IOLs can correct your astigmatism, and help you focus on things that are far away.



With Astigmatism



Without Astigmatism

These images are for illustrative purposes only and do not represent actual data derived from studies. These illustrative simulations are intended to help you better understand your vision in certain eye conditions. See more at [tecnisvisionsimulator.com](https://www.tecnisvisionsimulator.com).



Nearly 40% of people with cataracts also require astigmatism correction¹



When paired with a toric IOL, cataract surgery is a unique opportunity to address astigmatism, which typically occurs with other conditions

ASTIGMATISM



Distorts and blurs vision



Can lead to eye discomfort and headaches



Can be addressed with a toric IOL



Other solutions include prescription glasses and contacts made for astigmatism or an additional surgery



If you have been diagnosed with astigmatism, talk to your doctor about which lens option from TECNIS® IOLs might be right for your life.