



Increase patient and practice success by following the Fit Guide

ACUVUE® MULTIFOCAL PORTFOLIO WITH PUPIL OPTIMISED DESIGN TECHNOLOGY:



Offers a more **personalised solution** for your patients.**1



Provides a more precise fit to help keep the lens' optical design in the **right place**.¹

LENS DETAILS	1-DAY ACUVUE® MOIST MULTIFOCAL	ACUVUE® OASYS MULTIFOCAL 2-WEEKLY	ACUVUE® OASYS MAX 1-Day MULTIFOCAL
Material	etafilcon A	senofilcon A	senofilcon A
Diameter	14.3 mm	14.3 mm	14.3 mm
Base curve	8.4 mm	8.4 mm	8.4 mm
Technology	Embedded PVP*/LACREON® Technology	Embedded PVP*/HYDRACLEAR® PLUS Technology	TearStable™ Technology OptiBlue™ Light Filter‡
UV blocker*	Class 2	Class 1	Class 1
Dk/t*	25.5 x 10 ⁻⁹ §	147 x 10 ⁻⁹ §	147 x 10 ⁻⁹ §
Visibility tint	Yes	Yes	Yes Blue-green§2
Sphere	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW: +0.75D to +1.25D MID: +1.50D to +1.75D HIGH: +2.00D to +2.50D

Oxygen transmissibility at centre of a -3.00D lens using boundary-corrected, edge-corrected Dk values. Units: (cm/sec) (ml O2/ml x mm Hg) at 35°C. Dk determined via polarographic method.

+PVP=polyvinylpyrrolidone.



Visit the **ACUVUE® Multifocal Fitting Calculator** for quick & easy contact lens fitting & lens selection



** Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power. * All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00D lens. ‡ Filtering of High Energy Visible (HEV) light by contact lenses has not been demonstrated to confer any systemic and/or ocular health benefit to the user. The Eye Care Professional should be consulted for more information. § ACUVUE® OASYS MAX 1 Day has a unique blue-green appearance as a result of the combination of the blue-violet/high energy visible [HEV] light filter and the blue handling tint.

1. JJV Data on File 2022. ACUVUE® PUPIL OPTIMISED DESIGN TECHNOLOGY: JJVC Contact Lenses, Design Features, and Associated Benefits.
2. JJV Data on File 2022. TearStable™ Technology Definition.

ACUVUE® Contact Lenses are indicated for vision correction. For detailed product description and safety information, please consult the Instructions for Use or visit our Johnson & Johnson website www.e-ifu.com

© Johnson & Johnson Medical Ltd 2023. ACUVUE®, ACUVUE® MOIST, ACUVUE® OASYS, ACUVUE® OASYS MAX 1-Day MULTIFOCAL, LACREON®, HYDRACLEAR®, TearStable™, OptiBlue™ are registered trademarks of Johnson & Johnson. PP2023MLT7079

ACUVUE®
MULTIFOCAL
WITH PUPIL OPTIMISED DESIGN

ACUVUE®

YOUR QUICK AND EASY GUIDE TO FIT SUCCESS



Product images for illustrative purposes only



Unique
PUPIL OPTIMISED
DESIGN¹

ACUVUE®
MULTIFOCAL
Fit guide

Fit success
& patient
satisfaction

DESIGNED FOR SUPERIOR VISUAL PERFORMANCE.*1

Now available as both Daily Disposable and Reusable contact lenses.


ACUVUE®
MULTIFOCAL
WITH PUPIL OPTIMISED DESIGN

* Compared to prior JJV multifocal design; technology optimised for both the parameters of refractive error and add power for a multitude of viewing distances and light levels



The only brand with 100% of parameters
optimised by age & refraction1**

PUPIL OPTIMISED DESIGN



For illustrative purposes only. Pupil area can vary by ~20% at a given luminance.²

IN-BUILT PRECISION

ACUVUE® MULTIFOCAL PORTFOLIO with PUPIL OPTIMISED DESIGN provides a more **PRECISE FIT**: Hybrid Back Curve Technology **better matches the shape of the natural eye** to help keep the lens' optical design in the **right place**.¹



Product images for illustrative purposes only

Every parameter is designed to match different pupil sizes and provide the best balance of vision for that age and refraction¹

** Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power.
1 Across the power range of +6.00D to -9.00D.

INITIAL LENS SELECTION

1

Determine the Best Vision Sphere (BVS)

In the trial frame, confirm the least minus spherical prescription that provides the best distance VA[†]

2

Determine the sensory dominant eye

The +1.00D blur test recommended with the BVS in the trial frame rather than sighting methods.

3

Determine the multifocal ADD based upon the patient's needs

With the BVS in the trial frame, now determine the lowest ADD required to achieve good near vision.

Top Tip: Start with 0.50D less than the spectacle ADD and if necessary, increase in 0.25D steps until required near vision is achieved.

4

Select lens based on following tables

Multifocal ADD	Initial Lens Selection		Enhance Distance		Enhance Near	
	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye
+0.75 to +1.25	LOW	LOW	Use a spherical ACUVUE® lens	LOW	LOW	LOW & give extra +0.25D to dist. Rx
+1.50 to +1.75	MID	MID	LOW	MID	MID	MID & give extra +0.25D to dist. Rx
+2.00 to +2.50	MID	HIGH	MID	MID & give extra +0.25D to dist. Rx	MID	HIGH & give extra +0.25D to dist. Rx

AFTERCARE APPOINTMENT

If a **previously successful** ACUVUE® Multifocal wearer returns for an aftercare appointment reporting changes to distance or near vision:

- Do not over-refract their current lenses
- Do not use the 'Enhance distance / near' table above
- Do not adjust the ADD
- Remove lenses and repeat 'Initial lens selection' steps 1-4 above

[†] Proceed if astigmatism is less than or equal to 0.75DC. ⁹ Apply vertex distance correction if greater than +/- 4.00D.