



## Increase patient and practice success by following the Fit Guide

### ACUVUE® MULTIFOCAL PORTFOLIO WITH PUPIL OPTIMISED DESIGN TECHNOLOGY:



Offers a more **personalised solution** for your patients.\*\*1



Provides a more **precise fit** to help keep your patient's optics in the right place and the right shape<sup>1</sup>

LENS DETAILS	1-DAY ACUVUE® MOIST MULTIFOCAL	ACUVUE® OASYS MULTIFOCAL 2-WEEKLY	ACUVUE® OASYS MAX 1-Day MULTIFOCAL
Material	etafilcon A	senofilcon A	senofilcon A
Diameter	14.3 mm	14.3 mm	14.3 mm
Base curve	8.4 mm	8.4 mm	8.4 mm
Technology	Embedded PVP*/LACREON® Technology	Embedded PVP*/HYDRACLEAR® PLUS Technology	TearStable™ Technology OptiBlue™ Light Filter <sup>‡</sup>
UV blocker*	Class 2	Class 1	Class 1
Dk/t*	25.5 x 10 <sup>-9</sup>	147 x 10 <sup>-9</sup>	147 x 10 <sup>-9</sup>
Visibility tint	Yes	Yes	Yes   Blue-green <sup>§2</sup>
Sphere	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps)
ADD	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW +0.75D to +1.25D MID +1.50D to +1.75D HIGH +2.00D to +2.50D	LOW: +0.75D to +1.25D MID: +1.50D to +1.75D HIGH: +2.00D to +2.50D

# Oxygen transmissibility at centre of a -3.00D lens using boundary-corrected, edge-corrected Dk values. Units: (cm/sec) (ml O<sub>2</sub>/ml x mm Hg) at 35°C. Dk determined via polarographic method.

+PVP=polyvinylpyrrolidone.



Visit the **ACUVUE® Multifocal Fitting Calculator** for quick & easy contact lens fitting & lens selection



\*\* Compared to competitor's designs; technology optimised for both the parameters of refractive error and add power. \* All ACUVUE® contact lenses have Class 1 or Class 2 UV-blocking to help provide protection against transmission of harmful UV radiation to the cornea and into the eye. UV-absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV-absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. UV transmission measured with -1.00D lens. ‡ Filtering of High Energy Visible (HEV) light by contact lenses has not been demonstrated to confer any systemic and/or ocular health benefit to the user. The Eye Care Professional should be consulted for more information. § ACUVUE® OASYS MAX 1 Day has a unique blue-green appearance as a result of the combination of the blue-violet/high energy visible [HEV] light filter and the blue handling tint.

1. JJV Data on File 2022. ACUVUE® PUPIL OPTIMISED DESIGN TECHNOLOGY: JJVC Contact Lenses, Design Features, and Associated Benefits.

2. JJV Data on File 2022. TearStable™ Technology Definition.

ACUVUE® Contact Lenses are indicated for vision correction. For detailed product description and safety information, please consult the Instructions for Use or visit our Johnson & Johnson website [www.jnjvisioncare.ae](http://www.jnjvisioncare.ae)

ACUVUE®, ACUVUE® MOIST, ACUVUE® OASYS, ACUVUE® OASYS MAX 1-Day MULTIFOCAL, LACREON®, HYDRACLEAR®, TearStable™, OptiBlue™ are registered trademarks of Johnson & Johnson. © Johnson & Johnson Middle East, Inc. 2024. 2024PP15728

**ACUVUE®**  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN

ACUVUE®

## YOUR QUICK AND EASY GUIDE TO FIT SUCCESS



Product images for illustrative purposes only



Unique  
PUPIL OPTIMISED  
DESIGN<sup>1</sup>

ACUVUE®  
MULTIFOCAL  
Fit guide

Fit success  
& patient  
satisfaction

## DESIGNED FOR SUPERIOR VISUAL PERFORMANCE.\*1

Now available as both Daily Disposable and Reusable contact lenses.


ACUVUE®  
**MULTIFOCAL**  
WITH PUPIL OPTIMISED DESIGN

\* Compared to prior JJV multifocal design; technology optimised for both the parameters of refractive error and add power for a multitude of viewing distances and light levels



The only brand with 100% of parameters  
**optimised by age & refraction\*\*1**

## PUPIL OPTIMISED DESIGN



For illustrative purposes only. Pupil area can vary by ~20% at a given luminance.<sup>2</sup>

## IN-BUILT PRECISION

ACUVUE® MULTIFOCAL PORTFOLIO with PUPIL OPTIMISED DESIGN provides a more **PRECISE FIT**: Hybrid Back Curve Technology **better matches the shape of the natural eye** to help keep the lens' optical design in the **right place**.<sup>1</sup>



Product images for illustrative purposes only

**Every parameter is designed to match different pupil sizes  
and provide the best balance of vision for that age and refraction<sup>1</sup>**

\*\* Compared to competitors designs; technology optimised for both the parameters of refractive error and add power.  
° Across the power range of +6.00D to -9.00D.

## INITIAL LENS SELECTION

1

### Determine the Best Vision Sphere (BVS)

In the trial frame, confirm the least minus spherical prescription that provides the best distance VA<sup>†</sup>

2

### Determine the sensory dominant eye

The +1.00D blur test recommended with the BVS in the trial frame rather than sighting methods.

3

### Determine the lowest ADD based upon the patient's needs

With the BVS in the trial frame, now determine the lowest ADD required to achieve good near vision.

**Top Tip:** Start with 0.50D less than the spectacle ADD and if necessary, increase in 0.25D steps until required near vision is achieved.

4

### Select lens based on following tables

Spectacle ADD	Initial Lens Selection		Enhance Distance		Enhance Near	
	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye	Dominant Eye	Non-Dominant Eye
+0.75 to +1.25	LOW	LOW	Use a spherical ACUVUE® lens	LOW	LOW	LOW & give extra +0.25D to dist. Rx
+1.50 to +1.75	MID	MID	LOW	MID	MID	MID & give extra +0.25D to dist. Rx
+2.00 to +2.50	MID	HIGH	MID	MID & give extra +0.25D to dist. Rx	MID	HIGH & give extra +0.25D to dist. Rx

## AFTERCARE APPOINTMENT

If a **previously successful** ACUVUE® Multifocal wearer returns for an aftercare appointment reporting changes to distance or near vision:

- Do not over-refract their current lenses
- Do not use the 'Enhance distance / near' table above
- Do not adjust the ADD
- Remove lenses and repeat 'Initial lens selection' steps 1-4 above

<sup>†</sup> Proceed if astigmatism is less than or equal to 0.75DC. <sup>°</sup> Apply vertex distance correction if greater than +/- 4.00D.