EN- Instruction for Use ACUVUE Abiliti 1-Day Soft Therapeutic Lenses for Myopia Management

PLEASE READ CAREFULLY AND RETAIN FOR FUTURE REFERENCE. ACUVUE CONTACT LENSES (CL) ARE MEDICAL DEVICES AND SHOULD BE FITTED BY AN EYE CARE PROFESSIONAL (ECP). ALWAYS FOLLOW YOUR ECP'S DIRECTIONS AND THE INSTRUCTIONS IN THIS LEAFLET. IF YOU ARE THE PARENT OR GUARDIAN OF A CHILD WHO HAS BEEN PRESCRIBED TO WEAR, PLEASE EXPLAIN THE CONTENTS OF THIS INSTRUCTIONS FOR USE TO THE PATIENT AND ENSURE THAT THEY UNDERSTAND THIS INFORMATION.

1. PRODUCT DESCRIPTION and INTENDED USE: ACUVUE Abiliti 1-Day Soft Therapeutic Lenses for Myopia Management are intended to be worn for less than 24 hours while awake. Your ECP should prescribe the lenses and determine your wearing schedule. Your lenses do not require cleaning or disinfection and should be discarded upon removal. These Lenses are intended for Daily Wear for the optical correction of myopia (short-sightedness) while the lenses are worn and may slow the patient's prescription change over time. The contact lenses are intended to be used by children who, at the initiation of treatment, are 7 to 12 years old, with myopia between -0.75D and -4.50D and 1.00D or less astigmatism, with healthy eyes and are symptomatic regarding a need for distance vision correction.

Material: senofilcon A

Packaging Solution: Borate buffered saline with methyl ether cellulose. Borates (boric acid & sodium borate) are defined as CMR 1B substances in a concentration above 0.1% weight by weight and are safe when product is used according to label instructions.

WARNING: UV absorbing CL are not substitutes for UV protection. You should continue to use UV absorbing eyewear as directed.

CONTRAINDICATIONS: When wearing CL for vision correction, DO NOT USE CL in case of eye disease, injury, inflammation, infection, or abnormality.

• Any condition or disease of the eye or surrounding tissues induced or aggravated by wearing CL or use of CL solutions

• If eyes become red or irritated. Some people may develop an allergic response to chemicals or preservatives in contact lens care solutions.

WARNINGS: Failure to adhere to these warnings can lead to severe infection, vision loss or blindness. EYE PROBLEMS (INCLUDING CORNEAL ULCERS) CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION. IF YOU EXPERIENCE ANY OF THE FOLLOWING SYMPTOMS (Eye Discomfort, Loss of Vision, Vision Changes, Excessive Tearing, Eye Redness), you should IMMEDIATELY REMOVE YOUR CL AND PROMPTLY CONTACT YOUR ECP.

• Do not wear CL if there is any eye disease

• DO NOT WEAR YOUR CL WHEN SLEEPING

• DO NOT EXPOSE CL TO WATER, AS THIS MAY LEAD TO SERIOUS INFECTION. Discard and replace any CL exposed to water. (Do not wear CL during water sports).

• NEVER rinse your lenses in water from the tap.

Vision Quality and Potential Symptoms: Due to the intended optical design of these lenses, under certain circumstances, some patients may experience visual symptoms. These lenses will provide vision correction while worn, but the vision quality may not be as clear as with conventional soft contact lenses. Patients should exercise extra care if performing potentially hazardous, vision demanding activities. If the patient is not comfortable with how well they can see with these lenses, the patient should remove the lenses and contact the ECP.

PRECAUTIONS: DO NOT USE YOUR LENSES IF THE STERILE BLISTER PACKAGE IS OPENED, DAMAGED OR AFTER THE EXPIRY DATE.

• Only wear your own prescribed CL. Do not change your lens type or parameters without consulting your ECP. Ask your ECP about the follow-up visit schedule.

• ALWAYS CONTACT ECP BEFORE USING ANY MEDICINE IN YOUR EYES.

• NEVER use tweezers or other tools to remove your lenses from the lens container.

• Smoking and/or taking some medicines (like contraceptives or other medicines) may increase eye problems or eye discomfort.

• AVOID all harmful or irritating vapors/fumes. Keep your eyes closed during usage of aerosol products.

ADVERSE REACTIONS: Burning, stinging, itchy, and/or dry eyes - Reduced lens comfort or feeling of something in your eye - Swelling or inflammation in or around the eyes - Eye redness - Eyelid problems - Watery eyes and/ or unusual eye secretions - Poor or blurred vision - Rainbows or halos around objects- Sensitivity to light - Temporary impairment due to peripheral infiltrates, peripheral corneal ulcers, and corneal erosion. If you notice any problems, you should IMMEDIATELY REMOVE YOUR CL and contact your ECP.

DIRECTIONS FOR USE: Not following directions for use may result in severe eye problems.

• ALWAYS WASH YOUR HANDS BEFORE TOUCHING YOUR CL as CL are very sensitive to contamination.

• NEVER PUT YOUR LENSES IN YOUR MOUTH OR PUT SALIVA ON YOUR LENSES.

• NEVER ALLOW ANYONE ELSE TO WEAR YOUR LENSES.

• Follow your ECP instructions regarding wear and removal of your CL.

Lens Package: Carefully peel back the foil closure to reveal the lens.

Lens Insertion: Ensure the lens is not inside out. Place on your index finger to verify the shape. The lens should assume a natural, curved, bowl-like shape (Fig. 1). With the lens on your index finger, use your other hand to hold your upper eyelid so you will not blink. Place the lens onto the eye. (Fig 2 &3).

Lens Removal: Look up, slide the lens to the lower part of the eye using the forefinger. Gently pinch the lens between the thumb and forefinger to remove the lens. (Fig. 4)

EMERGENCIES: If chemicals of any kind (household products, gardening solutions, laboratory chemicals, etc.) are splashed into your eyes: FLUSH EYES IMMEDIATELY WITH TAP WATER AND IMMEDIATELY CONTACT YOUR ECP OR VISIT A HOSPITAL EMERGENCY ROOM RIGHT AWAY.

REPORTING OF ADVERSE REACTIONS (Side Effects): Any incident experienced whilst wearing ACUVUE Brand CL should be reported to the manufacturer and/or its authorized representative and/or to your national authority.

EN- The following symbols may appear on the labels or packaging

Symbol	Definition
	Date of Manufacture
	Manufacture
	Use-by Date (expiration date)
	Batch Code
	Do Not Use if Package is Damaged
	Consult Instructions for Use or Consult Electronic Instructions for use
	Caution
	Sterilized Using Steam Heat
	Do Not Re-Use (Single Use)
	Medical Device in the European Community
	Indicates a Single Sterile Barrier System
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner
	UV Blocking
	Authorized Representative in the European Community
	CE-mark and Identification Number of Notified Body
	Diameter
	Base Curve
	Diopter (lens power)
	Contains Hazardous Substances

FR - Guide d'utilisation des lentilles ACUVUE Abiliti 1-Day Soft Therapeutic Lenses

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©Johnson & Johnson Vision Care Companies
USA: Johnson & Johnson Vision Care, Inc., 7500 Centurion Parkway, Jacksonville, Florida, 32256, USA