


TECNIS Simplicity™ — The Next-Generation TECNIS® IOL Disposable Delivery System

Simplifying delivery of the IOL you prefer.



BENEFITS:



Easy 3-step process -
hydrate, advance, deliver¹



Advanced design for a
smooth and controlled
IOL delivery¹



Flexibility to use
Balanced Salt Solution
or OVD for hydration²



For more information on the TECNIS Simplicity™ Delivery System,
please contact your Johnson & Johnson Vision representative.

TECNIS SIMPLICITY™
DELIVERY SYSTEM

Johnson & Johnson VISION

The Next-Generation TECNIS® IOL Disposable Delivery System



**Model:
DIB00**



**Model:
DIU100-800**

OPTIC CHARACTERISTICS ^{2,3}	
Power Range:	+5.0 D to +34.0 D in 0.5 diopter increments
Model Numbers:	DIB00
Cylinder Powers - IOL Plane	n/a
Cylinder Powers - Corneal Plane	n/a
Diameter:	6.0 mm
Shape:	Biconvex, anterior aspheric surface
Material:	UV-absorbing hydrophobic acrylic
Refractive Index:	1.47 at 35°C
Edge Design:	ProTEC 360° Edge; frosted edge design; continuous posterior square edge
OPTICAL BIOMETRY/APPLANATION ULTRASOUND BIOMETRY	
A-constant:	119.3*
Theoretical AC Depth: ⁴	5.7 mm*
Surgeon Factor: ⁴	1.96 mm*
HAPTIC CHARACTERISTICS ^{2,3}	
Lens Overall Diameter:	13.0 mm
Style:	C
Design:	Tri-Fix, haptics offset from optics, 1-piece lens
Incision size:	2.2 mm – 2.4 mm

*Derived from clinical evaluation results of the TECNIS® 1-Piece IOL Platform for optical biometry.

*A-Constant theoretically derived for ultrasound biometry.

+5.0 D to +34.0 D in 0.5 diopter increments									
DIU100	DIU150	DIU225	DIU300	DIU375	DIU450	DIU525	DIU600	DIU700	DIU800
1.00 D	1.50 D	2.25 D	3.00 D	3.75 D	4.50 D	5.25 D	6.00 D	7.00 D	8.00
0.69 D	1.03 D	1.54 D	2.06 D	2.57 D	3.08 D	3.60 D	4.11 D	4.80 D	5.48 D
6.0 mm									
Biconvex, anterior aspheric surface									
UV-absorbing hydrophobic acrylic									
1.47 at 35°C									
ProTEC 360° Edge; frosted edge design; continuous posterior square edge									
118.8†									
5.4 mm									
1.68 mm									
13.0 mm									
C									
Tri-Fix, haptics offset from optics, 1-piece lens									
2.2 mm – 2.4 mm									

References: 1. Delta IOL Delivery System - Human Factors Validation Report. Dec. 12, 2018. REF2019CT4449. 2. TECNIS Eyhance™ IOL with TECNIS Simplicity™ Delivery System, Model DIB00 - DfU INT - Z311246P, Rev. B, 08/2019. REF2019OTH4409. 3. TECNIS Eyhance™ Toric II IOL with TECNIS Simplicity™ Delivery System, Model DIU - DfU INT - Z311515P, Rev. A, 09/2020. REF2020CT4293. 4. Holladay JT. International Intraocular Lens & Implant Registry 2004. *J Cataract Refract Surg* 2004;30(1):207-229. REF2019CT4464.

For healthcare professionals only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.