



Simplifying  
IOL Delivery.  
Improving  
Efficiency.

Johnson & Johnson

**TECNIS** SIMPLICITY™  
DELIVERY SYSTEM

TECNIS SIMPLICITY™  
DELIVERY SYSTEM

# The Safety of a Single-Use System. For the Delivery of High Quality Lenses.

The preloaded, fully disposable **TECNIS Simplicity™** Delivery System is designed to help to prevent manual loading errors, streamlines lens delivery, and protects against contamination to help elevate the standard of care.<sup>1</sup>


**TECNIS Simplicity™** comes conveniently preloaded with select **TECNIS™** IOLs.





## Improved Efficiency

Simple 3-step process:  
hydrate, advance  
and deliver<sup>2</sup>



## Smooth and Controlled Delivery

Designed to safely provide  
placement of the IOL into  
the capsular bag<sup>2</sup>



## Flexibility

Hydration can be done  
using Balanced Salt  
Solution or OVD<sup>1</sup>

The promoted product is a medical device. For Healthcare Professionals Only. Please reference the Instructions for Use for a complete list of Indications and Important Safety Information and contact our specialists in case of any question.

REFERENCES: 1. TECNIS™ 1-Piece IOL with TECNIS Simplicity™ Delivery System, Model DCB00, elfU, Z311526E, current revision. 2. Delta IOL Delivery System - Human Factors Validation Report. Dec. 12, 2018. REF2019CT4449.

**TECNIS** SIMPLICITY™  
DELIVERY SYSTEM

# Simplifying IOL Delivery.

Choose the next-generation TECNIS™ IOL delivery system for:



## Improved Efficiency

Simple 3-step process:  
hydrate, advance and deliver<sup>2</sup>



## Eliminated Manual Loading

Helps to prevent manual loading errors and no IOL touches, reducing haptic breakage, marking and crimping<sup>1</sup>



## Enhanced Safety

Minimized risk of infection associated with contamination<sup>1</sup>



## Flexibility

Hydration can be done using  
Balanced Salt Solution or OVD<sup>1</sup>



## Smooth and Controlled Delivery


Designed to safely provide placement of the IOL into the capsular bag<sup>2</sup>



## Convenience

Preloaded and fully disposable<sup>1</sup>

# Enhanced Ergonomics and Optimized Features Designed to Deliver on Efficiency.



## Enlarged Button

More surgeon control<sup>2</sup>

## Ergonomic Design

Sturdy, comfortable grip; better grasp<sup>2</sup>

## Firmer Plunger

Smooth and controlled delivery<sup>2</sup>

## Lens Case

Preloaded with select TECNIS™ 1-piece IOLs<sup>1</sup>

## Hydration Port

For Balanced Salt Solution or OVD<sup>1</sup>


## Beveled Tip

For a 2.2–2.4 mm incision<sup>1</sup>

# The *Simplicity* of a 3-Step Process.



**STEP 1: HYDRATE**



**STEP 2: ADVANCE**



**STEP 3: DELIVER**

**With the Flexibility to Hydrate How You Prefer.**

TECNIS Simplicity™ is a preloaded IOL delivery system that allows you to choose from 2 hydration methods:

- ✓ Balanced Salt Solution
- ✓ OVD