















TECHNICAL SPECIFICATIONS

TECHNICAL SPECIFICATIONS	 ACUVUE® OASYS MAX 1-Day	 ACUVUE® OASYS MAX 1-Day MULTIFOCAL	 ACUVUE OASYS® 1-DAY	 ACUVUE OASYS® 1-DAY for ASTIGMATISM	 1-DAY ACUVUE® MOIST	 1-DAY ACUVUE® MOIST for ASTIGMATISM	 1-DAY ACUVUE® MOIST MULTIFOCAL	 1-DAY ACUVUE® DEFINE®	 1-DAY ACUVUE® DEFINE® FRESH	 ACUVUE® OASYS®	 ACUVUE® OASYS® for ASTIGMATISM	 ACUVUE® OASYS MULTIFOCAL	 ACUVUE® 2
Base curve	8.5mm, 9.0mm	8.4mm	8.5mm, 9.0mm	8.5mm	8.5mm, 9.0mm	8.5mm	8.4mm	8.5mm	8.5mm	8.4mm	8.6mm	8.4mm	8.3mm, 8.7mm
Diameter	14.3mm	14.3mm	14.3mm	14.3mm	14.2mm	14.5mm	14.3mm	14.2mm	14.2mm	14.0mm	14.5mm	14.3mm	14.0mm
 Approximate UV Blocking**	100% UVB, 99.9% UVA	100% UVB, 99.9% UVA	Blocks-99.9% of UVB, 96% of UVA	Blocks-99.9% of UVB, 96% of UVA	Blocks-99% of UVB, 85% of UVA	Blocks-99% of UVB, 85% of UVA	Blocks-99% of UVB, 86% of UVA	Blocks-99% of UVB, 85% of UVA	Blocks-99% of UVB, 85% of UVA	Blocks-99.9% of UVB, 96% of UVA	Blocks-99.9% of UVB, 96% of UVA	Blocks-99.9% of UVB, 96.1% of UVA	Blocks-99% of UVB, 85 % of UVA
Recommended Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear 2-Weekly Replacement
Power Range (D)	+0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps) -0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps)	-9.00D to +6.00D (0.25D steps) Low +0.75D to +1.25D Mid +1.50D to +1.75D High +2.00D to +2.50D	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Low minus 0.00D to -6.00D (0.25D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10° - 180°) CYL: -2.25 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°) High minus -6.50D to -9.00D (0.50D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160, 170) Plus +0.25D to +4.00D (0.25D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)	Low minus 0.00 to -6.00D (0.25D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10° - 180°) CYL: -2.25 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°) High minus -6.50D to -9.00D (0.50D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160, 170) Plus +0.25D to +4.00D (0.25D steps) CYL: -0.75, -1.25, -1.75 (AXIS: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°)	Sphere -9.00D to +6.00D (0.25D steps) ADDs Low (+0.75D to +1.25D) Mid (+1.50D to +1.75D) High (+2.00D to +2.50D)	-0.50 to -6.00D (0.25 steps) -6.50 to -9.00D (0.50 steps) 0.00 in VIVID STYLE only	-0.50D to -6.00D (0.25 steps) -6.50D to -9.00D (0.50 steps)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Low minus: 0.00 to -6.00D (0.25D steps) CYL: -0.75, -1.25, -1.75, -2.25, -2.75 AXIS: 10° to 180° (10° steps) High minus: -6.50D to -9.00D (0.50D steps) CYL: -0.75, -1.25, -1.75, -2.25, -2.75 AXIS: 10° to 180° (10° steps) Plus: +0.25D to +6.00D (0.25D steps) CYL: -0.75, -1.25, -1.75, -2.25, -2.75 AXIS: 10° to 180° (10° steps)	Sphere -9.00 to +6.00D (0.25D steps) ADDs Low (+0.75D to +1.25D) Mid (+1.50D to +1.75D) High (+2.00D to +2.50D)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)

ACUVUE®, ACUVUE® OASYS MAX 1-Day and ACUVUE® OASYS MAX 1-Day MULTIFOCAL, ACUVUE® OASYS® 1-DAY, ACUVUE® OASYS® 1-DAY for ASTIGMATISM, ACUVUE® OASYS MULTIFOCAL, ACUVUE® OASYS®, ACUVUE® OASYS® for ASTIGMATISM, 1-DAY ACUVUE® MOIST, 1-DAY ACUVUE® MOIST for ASTIGMATISM, 1-DAY ACUVUE® MOIST MULTIFOCAL, 1-DAY ACUVUE® DEFINE, ACUVUE® 2 are trademarks of Johnson & Johnson Middle East Inc, © Johnson & Johnson Middle East Inc.2025.  
For more information on proper wear, care and safety, please consult the Instructions for Use or visit our J&J website [www.jnjvisioncare.ae](http://www.jnjvisioncare.ae)  
2025PP06757