

TECHNICAL SPECIFICATIONS

ACUVUE®, ACUVUE® OASYS MAX 1-Day, ACUVUE® OASYS MAX 1-Day MULTIFOCAL, ACUVUE® OASYS MAX 1-Day for ASTIGMATISM, ACUVUE® OASYS® 1-DAY, ACUVUE® OASYS® 1-DAY for ASTIGMATISM, ACUVUE® OASYS MULTIFOCAL, ACUVUE® OASYS®, ACUVUE® OASYS® for ASTIGMATISM, 1-DAY ACUVUE® MOIST, 1-DAY ACUVUE® MOIST for ASTIGMATISM, 1-DAY ACUVUE® MOIST MULTIFOCAL and 1-DAY ACUVUE® DEFINE, ACUVUE® 2 are trademarks of Johnson & Johnson Middle East Inc, © Johnson & Johnson Middle East Inc. 2026. For more information on proper wear, care and safety, please consult the Instructions for Use or visit our Johnson & Johnson Vision website www.jnjvisionpro.com/en-me/ 2025PP13679.

Johnson&Johnson

ACUVUE®

EYE-INSPIRED™

INNOVATIONS¹

Every ACUVUE® contact lens is inspired by the natural function of the eye¹ based on a deep understanding of patient needs.



MOISTURE-INFUSED¹
FOR HYDRATION



INVISIBLE EDGE¹
FOR COMFORT



UV BLOCKING¹
FOR EYE HEALTH



ANTIOXIDANTS¹
TO HELP SUPPORT COMFORT



EYELID STABILISED¹
FOR ASTIGMATISM



PUPIL OPTIMISED¹
FOR MULTIFOCAL VISION

In ACUVUE® Contact Lenses ^

In ACUVUE® contact lenses for Astigmatism

In ACUVUE® Multifocal contact lenses

TECHNICAL SPECIFICATIONS															
	ACUVUE® OASYS MAX 1-Day	ACUVUE® OASYS MAX 1-Day ASTIGMATISM	ACUVUE® OASYS MAX 1-Day MULTIFOCAL	ACUVUE OASYS® 1-DAY	ACUVUE OASYS® 1-DAY for ASTIGMATISM	1-DAY ACUVUE® MOIST	1-DAY ACUVUE® MOIST for ASTIGMATISM	1-DAY ACUVUE® MOIST MULTIFOCAL	1-DAY ACUVUE® DEFINE™	1-DAY ACUVUE® DEFINE™ FRESH	ACUVUE® OASYS®	ACUVUE® OASYS® for ASTIGMATISM	ACUVUE® OASYS MULTIFOCAL	ACUVUE® 2	
	SPHERE	TORIC	MULTIFOCAL	SPHERE	TORIC	SPHERE	TORIC	MULTIFOCAL	SPHERE	SPHERE	SPHERE	TORIC	MULTIFOCAL	SPHERE	
Base curve	8.5mm, 9.0mm	8.5mm	8.4mm	8.5mm, 9.0mm	8.5mm	8.5mm, 9.0mm	8.5mm	8.4mm	8.5mm	8.5mm	8.4mm, 8.8mm	8.6mm	8.4mm	8.3mm, 8.7mm	
Diameter	14.3mm	14.3mm	14.3mm	14.3mm	14.3mm	14.2mm	14.5mm	14.3mm	14.2mm	14.2mm	14.0mm	14.5mm	14.3mm	14.0mm	
 Approximate UV Blocking**	100% UVB, 99.9% UVA	100% UVB, 99.9% UVA	100% UVB, 99.9% UVA	99.9% of UVB, 96% of UVA	99.9% of UVB, 96% of UVA	99% of UVB, 85% of UVA	99% of UVB, 85% of UVA	99% of UVB, 86% of UVA	99% of UVB, 85% of UVA	99% of UVB, 85% of UVA	99.9% of UVB, 96.1% of UVA	99.9% of UVB, 96.1% of UVA	99.9% of UVB, 96.1% of UVA	99% of UVB, 88 % of UVA	
Recommended Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Single use 1-Day Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear: 2-Weekly Replacement Extended wear: 1-Weekly Replacement	Daily wear: 2-Weekly Replacement	
Power Range (D)	0.00 (Plano) -0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Low minus 0.00 to -6.00D (0.25D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10° to 180° (in 10° steps) Cyl: -2.25 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° High minus -6.50D to -9.00D (0.50D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° Plus +0.25D to +4.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps) ADDS: LOW: +0.75 to +1.25 MID: +1.50 to +1.75 HIGH: +2.00 to +2.50	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Low minus 0.00D to -6.00D (0.25D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10° to 180° (in 10° steps) Cyl: -2.25 Axis: , 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° High minus -6.50D to -9.00D (0.50D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160, 170, 180 Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° Plus +0.25D to +4.00D (0.25D steps)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)	Low minus 0.00 to -6.00D (0.25D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10° to 180° (in 10° steps) Cyl: -2.25 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° High minus -6.50D to -9.00D (0.50D steps) Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160, 170, 180 Cyl: -0.75, -1.25, -1.75 Axis: 10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180° Plus +0.25D to +4.00D (0.25D steps)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps)	-9.00D to +6.00D (0.25D steps) ADDS: LOW: +0.75 to +1.25 MID: +1.50 to +1.75 HIGH: +2.00 to +2.50	-0.50 to -6.00D (0.25D steps) -6.50 to -9.00D (0.50D steps) 0.00 in VIVID STYLE only	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)	Low minus 0.00 to -6.00D (0.25D steps) Cyl: -0.75, -1.25, -1.75, -2.25, -2.75 Axis: 10° to 180° (in 10° steps) High minus -6.50D to -9.00D (0.50D steps) Cyl: -0.75, -1.25, -1.75, -2.25, -2.75 Axis: 10° to 180° (in 10° steps) Cyl: -0.75, -1.25, -1.75, -2.25, -2.75 Axis: 10° to 180° (in 10° steps) Plus +0.25D to +6.00D (0.25D steps)	-9.00 to +6.00D (0.25D steps) ADDS LOW: +0.75 to +1.25 MID: +1.50 to +1.75 HIGH: +2.00 to +2.50	-0.50D to -6.00D (0.25D steps) -6.50D to -12.00D (0.50D steps) +0.50D to +6.00D (0.25D steps) +6.50D to +8.00D (0.50D steps)